

Antibodies Product Datasheet

Product Name:	ERK5 Mouse mAb(6D10)
Catalog Number:	N1503
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ERK5
Concentration:	1 mg/ml
Gene Name:	MAPK7
Protein Name:	ERK5
Other Name:	
Human Swiss Prot No.:	Q13164
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.79% sodium azide.
Specificity:	This antibody detects endogenous levels of ERK5 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year