

Antibodies Product Datasheet

Product Name:	PCSK9 Mouse mAb(5G4)
Catalog Number:	N1505
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human PCSK9
Concentration:	1 mg/ml
Gene Name:	PCSK9 NARC1 PSEC0052
Protein Name:	PCSK9
Other Name:	Proprotein convertase subtilisin/kexin type 9 (EC 3.4.21.-) (Neural apoptosis-regulated convertase 1) (NARC-1) (Proprotein convertase 9) (PC9) (Subtilisin/kexin-like protease PC9)
Human Swiss Prot No.:	Q8NBP7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.84% sodium azide.
Specificity:	This antibody detects endogenous levels of PCSK9 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year