

Antibodies Product Datasheet

Product Name:	GAP-43 Monoclonal Antibody(Mix)
Catalog Number:	N1512
Reactivity:	Human,Mouse,Rat
Applications:	WB,IHC-p,IF(paraffin section)
Dilutions:	WB: 1:1000-2000 IHC:1:200-500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of GAP-43
Concentration:	
Gene Name:	GAP43
Protein Name:	Neuromodulin
Other Name:	GAP43; Neuromodulin; Axonal membrane protein GAP-43; Growth-associated protein 43; Neural phosphoprotein B-50; pp46
Human Swiss Prot No.:	P17677
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous GAP-43 protein.
Storage and Stability:	-20°C/1 year