

Antibodies Product Datasheet

Product Name:	MICU1 Monoclonal Antibody(Mix)
Catalog Number:	N1518
Reactivity:	Human,Mouse,Rat
Applications:	WB,IF/ICC,IHC-p
Dilutions:	WB: 1:1000-2000 IHC: 1:100-200 IF 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of MICU1
Concentration:	
Gene Name:	MICU1
Protein Name:	Calcium uptake protein 1, mitochondrial
Other Name:	Calcium uptake protein 1, mitochondrial; Atopy-related autoantigen CALC; ara CALC; Calcium-binding atopy-related autoantigen 1; allergen Hom s 4
Human Swiss Prot No.:	Q9BPX6
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous MICU1 protein.
Storage and Stability:	-20°C/1 year