

Antibodies Product Datasheet

Product Name:	Active Caspase-3 Monoclonal Antibody(5E1)
Catalog Number:	N1519
Reactivity:	Human,Mouse,Rat,Chicken
Applications:	IF/ICC,WB,IHC-p
Dilutions:	IF: 1:50-200 WB: 1:500-1000 IHC: 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of Active Caspase-3
Concentration:	
Gene Name:	CASP3
Protein Name:	Caspase3
Other Name:	CASP3; CPP32; Caspase-3; CASP-3; Apopain; Cysteine protease CPP32; CPP-32; Protein Yama; SREBP cleavage activity 1; SCA-1
Human Swiss Prot No.:	P42574
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Active Caspase-3 protein.
Storage and Stability:	-20°C/1 year