

Antibodies Product Datasheet

Product Name:	Bovine Serum Albumin Monoclonal Antibody(Mix-mA)
Catalog Number:	N1521
Reactivity:	Cow
Applications:	WB
Dilutions:	WB: 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified Protein
Concentration:	
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Bovine Serum Albumin protein.
Storage and Stability:	-20°C/1 year