
Antibodies Product Datasheet

Product Name:	Acetyl Lysine Monoclonal Antibody(10B10)
Catalog Number:	N1522
Reactivity:	Species independent
Applications:	WB,IHC-p,IF(paraffin section) IP
Dilutions:	WB: 1:1000-2000 IHC: 1:200-500 IP 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified Protein
Concentration:	
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Acetyl Lysine protein.
Storage and Stability:	-20°C/1 year