

Antibodies Product Datasheet

Product Name:	HSP90 α Monoclonal Antibody(1F6)
Catalog Number:	N1531
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF(paraffin section)
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the C-terminal region of human HSP90A. at AA range: 660-740
Concentration:	1 mg/ml
Gene Name:	HSP90AA1
Protein Name:	Heat shock protein HSP 90-alpha (Heat shock 86 kDa) (HSP 86) (HSP86) (Renal carcinoma antigen NY-REN-38)
Other Name:	Heat shock protein HSP 90-alpha (Heat shock 86 kDa) (HSP 86) (HSP86) (Renal carcinoma antigen NY-REN-38)
Human Swiss Prot No.:	P07900
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous HSP90 α protein
Storage and Stability:	-20°C/1 year