

Antibodies Product Datasheet

Product Name:	PI3 Kinase P85 α Monoclonal Antibody(3B7)
Catalog Number:	N1540
Reactivity:	Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of PI3 Kinase P85 α
Concentration:	1 mg/ml
Gene Name:	PIK3R1
Protein Name:	Phosphatidylinositol 3-kinase regulatory subunit alpha (PI3-kinase regulatory subunit alpha) (PI3K regulatory subunit alpha) (PtdIns-3-kinase regulatory subunit alpha) (Phosphatidylinositol 3-kinase 85 kDa regulatory subunit alpha) (PI3-kinase subunit p85
Other Name:	PIK3R1; GRB1; Phosphatidylinositol 3-kinase regulatory subunit alpha; PI3-kinase regulatory subunit alpha; PI3K regulatory subunit alpha; PtdIns-3-kinase regulatory subunit alpha; Phosphatidylinositol 3-kinase 85 kDa regulatory subunit alpha; PI3-kinase s
Human Swiss Prot No.:	P27986
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PI3 Kinase P85 α protein detects endogenous levels of PI3 Kinase P85 α
Storage and Stability:	-20°C/1 year