

Antibodies Product Datasheet

Product Name:	TBP/TATA Binding Protein Monoclonal Antibody(4H2)
Catalog Number:	N1548
Reactivity:	Human,Rat,Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of TBP/TATA Binding Protein
Concentration:	1 mg/ml
Gene Name:	TBP
Protein Name:	TATA-box-binding protein (TATA sequence-binding protein) (TATA-binding factor) (TATA-box factor) (Transcription initiation factor TFIID TBP subunit)
Other Name:	TBP; GTF2D1; TF2D; TFIID; TATA-box-binding protein; TATA sequence-binding protein; TATA-binding factor; TATA-box factor; Transcription initiation factor TFIID TBP subunit
Human Swiss Prot No.:	P20226
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TBP/TATA Binding Protein protein detects endogenous levels of TBP/TATA Binding Protein
Storage and Stability:	-20°C/1 year