

Antibodies Product Datasheet

Product Name:	Cyclophilin B Monoclonal Antibody(2B10)
Catalog Number:	N1549
Reactivity:	Human,Rat,Mouse
Applications:	IF/ICC, WB, IHC-p
Dilutions:	IF: 1:50-200 WB 1:1000-2000 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Cyclophilin B
Concentration:	1 mg/ml
Gene Name:	PPIB
Protein Name:	Peptidyl-prolyl cis-trans isomerase B (PPIase B) (EC 5.2.1.8) (CYP-S1) (Cyclophilin B) (Rotamase B) (S-cyclophilin) (SCYLP)
Other Name:	PPIB; CYPB; Peptidyl-prolyl cis-trans isomerase B; PPIase B; CYP-S1; Cyclophilin B; Rotamase B; S-cyclophilin; SCYLP
Human Swiss Prot No.:	P23284
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cyclophilin B protein detects endogenous levels of Cyclophilin B
Storage and Stability:	-20°C/1 year