

Antibodies Product Datasheet

Product Name:	HSPB8/HSP22 Monoclonal Antibody(2C3)
Catalog Number:	N1562
Reactivity:	Human,Rat,Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of HSPB8/HSP22
Concentration:	1 mg/ml
Gene Name:	HSPB8
Protein Name:	Heat shock protein beta-8 (HspB8) (Alpha-crystallin C chain) (E2-induced gene 1 protein) (Protein kinase H11) (Small stress protein-like protein HSP22)
Other Name:	HSPB8; CRYAC; E2IG1; HSP22; PP1629; Heat shock protein beta-8; HspB8; Alpha-crystallin C chain; E2-induced gene 1 protein; Protein kinase H11; Small stress protein-like protein HSP22
Human Swiss Prot No.:	Q9UJY1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HSPB8/HSP22 protein detects endogenous levels of HSPB8/HSP22
Storage and Stability:	-20°C/1 year