

Antibodies Product Datasheet

Product Name:	α -tubulin (Acetyl Lys40) Monoclonal Antibody(5H5)
Catalog Number:	N1575
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of α -tubulin (Acetyl Lys40)
Concentration:	1 mg/ml
Gene Name:	TUBA1B
Protein Name:	Tubulin alpha-1B chain (Alpha-tubulin ubiquitous) (Tubulin K-alpha-1) (Tubulin alpha-ubiquitous chain)
Other Name:	Tubulin alpha-1B chain (Alpha-tubulin ubiquitous) (Tubulin K-alpha-1) (Tubulin alpha-ubiquitous chain)
Human Swiss Prot No.:	P68363
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	α -tubulin (Acetyl Lys40) protein detects endogenous levels of α -tubulin (Acetyl Lys40)
Storage and Stability:	-20°C/1 year