

Antibodies Product Datasheet

Product Name:	Digoxin Monoclonal Antibody(4F6)
Catalog Number:	N1584
Reactivity:	
Applications:	ELISA
Dilutions:	ELISA 1/10000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Digoxin
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	
Specificity:	
Storage and Stability:	-20°C/1 year