

## Antibodies Product Datasheet

|                               |                                                                                                                                                                                                                        |
|-------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | PPAR Delta mouse Monoclonal Antibody(2F9)                                                                                                                                                                              |
| <b>Catalog Number:</b>        | N1585                                                                                                                                                                                                                  |
| <b>Reactivity:</b>            | Human,Rat,Mouse                                                                                                                                                                                                        |
| <b>Applications:</b>          | IF/ICC,IHC-p                                                                                                                                                                                                           |
| <b>Dilutions:</b>             | IF: 1:50-200 IHC 1:100-200                                                                                                                                                                                             |
| <b>Modification:</b>          | Unmodified/ Total                                                                                                                                                                                                      |
| <b>Source:</b>                | Monoclonal Mouse                                                                                                                                                                                                       |
| <b>Purification:</b>          | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.                                                                                                             |
| <b>Immunogen:</b>             | Recombinant Protein of PPAR Delta of PPAR Delta                                                                                                                                                                        |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                                                                                |
| <b>Gene Name:</b>             | PPARD NR1C2 PPARG                                                                                                                                                                                                      |
| <b>Protein Name:</b>          | PPAR Delta                                                                                                                                                                                                             |
| <b>Other Name:</b>            | Peroxisome proliferator-activated receptor delta (PPAR-delta) (NUC1) (Nuclear hormone receptor 1) (NUC1) (Nuclear receptor subfamily 1 group C member 2) (Peroxisome proliferator-activated receptor beta) (PPAR-beta) |
| <b>Human Swiss Prot No.:</b>  | Q03181                                                                                                                                                                                                                 |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                |
| <b>Specificity:</b>           | PPAR Delta protein detects endogenous levels of PPAR Delta                                                                                                                                                             |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                                                                           |