

Antibodies Product Datasheet

Product Name:	IL-8 mouse Monoclonal Antibody(13F8)
Catalog Number:	N1593
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic C-TERM Peptide of IL-8
Concentration:	1 mg/ml
Gene Name:	IL8
Protein Name:	IL8
Other Name:	IL8
Human Swiss Prot No.:	P10145
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-8 protein detects endogenous levels of IL8
Storage and Stability:	-20°C/1 year