

Antibodies Product Datasheet

Product Name:	a-actinin Monoclonal Antibody(4A12)
Catalog Number:	N1594
Reactivity:	Human,Rat,Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of a-actinin of ACTN1
Concentration:	1 mg/ml
Gene Name:	ACTN1
Protein Name:	ACTN1
Other Name:	ACTN1
Human Swiss Prot No.:	P12814
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	a-actinin Monoclonal Antibody(4A12) detects endogenous levels of ACTN1
Storage and Stability:	-20°C/1 year