

Antibodies Product Datasheet

Product Name:	Bax mouse Monoclonal Antibody(6F11)
Catalog Number:	N1601
Reactivity:	Human,Rat,Mouse,Chicken
Applications:	IF/ICC, WB, IHC-p
Dilutions:	IF: 1:50-200 WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Bax at AA range of 30-110
Concentration:	1 mg/ml
Gene Name:	BAX
Protein Name:	BAX
Other Name:	BAX
Human Swiss Prot No.:	Q07812
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Bax protein detects endogenous levels of BAX
Storage and Stability:	-20°C/1 year