

Antibodies Product Datasheet

Product Name:	LC3B mouse Monoclonal Antibody(9H5)
Catalog Number:	N1624
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of LC3B of MAP1LC3B
Concentration:	1 mg/ml
Gene Name:	MAP1LC3B
Protein Name:	MAP1LC3B
Other Name:	MAP1LC3B
Human Swiss Prot No.:	Q9GZQ8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LC3B protein detects endogenous levels of MAP1LC3B
Storage and Stability:	-20°C/1 year