

Antibodies Product Datasheet

Product Name:	Luciferase mouse Monoclonal Antibody(6B8)
Catalog Number:	N1628
Reactivity:	Firefly
Applications:	IF/ICC, WB,
Dilutions:	IF: 1:50-200 WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Luciferase of Luciferase
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	Luciferase
Other Name:	
Human Swiss Prot No.:	P08659
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Luciferase protein detects endogenous levels of Luciferase
Storage and Stability:	-20°C/1 year