

Antibodies Product Datasheet

Product Name:	Epsilon Tubulin mouse Monoclonal Antibody(2B3)
Catalog Number:	N1636
Reactivity:	Human,Rat,Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Epsilon Tubulin
Concentration:	1 mg/ml
Gene Name:	TUBE1 TUBE
Protein Name:	Epsilon Tubulin
Other Name:	Tubulin epsilon chain (Epsilon-tubulin)
Human Swiss Prot No.:	Q9UJT0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Epsilon Tubulin protein detects endogenous levels of Epsilon Tubulin
Storage and Stability:	-20°C/1 year