

## Antibodies Product Datasheet

|                               |                                                                                                            |
|-------------------------------|------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | CHOP mouse Monoclonal Antibody(2B1)                                                                        |
| <b>Catalog Number:</b>        | N1650                                                                                                      |
| <b>Reactivity:</b>            | Human,Rat,Mouse                                                                                            |
| <b>Applications:</b>          | WB,IF/ICC,IHC-p                                                                                            |
| <b>Dilutions:</b>             | WB 1:1000-2000, IHC 1:100-200 IF 1:200                                                                     |
| <b>Modification:</b>          | Unmodified/ Total                                                                                          |
| <b>Source:</b>                | Monoclonal Mouse                                                                                           |
| <b>Purification:</b>          | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen. |
| <b>Immunogen:</b>             | Synthetic Peptide of CHOP at AA range of 10-90                                                             |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                    |
| <b>Gene Name:</b>             | DDIT3                                                                                                      |
| <b>Protein Name:</b>          | DDIT3                                                                                                      |
| <b>Other Name:</b>            | DDIT3                                                                                                      |
| <b>Human Swiss Prot No.:</b>  | P35638                                                                                                     |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                    |
| <b>Specificity:</b>           | CHOP protein detects endogenous levels of DDIT3                                                            |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                               |