

Antibodies Product Datasheet

Product Name:	ERK1 mouse Monoclonal Antibody(5E9)
Catalog Number:	N1656
Reactivity:	Human,Rat,Mouse
Applications:	IF/ICC,IHC-p
Dilutions:	IF: 1:50-200 IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of ERK1 of MAPK3
Concentration:	1 mg/ml
Gene Name:	MAPK3
Protein Name:	MAPK3
Other Name:	MAPK3
Human Swiss Prot No.:	P27361
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ERK1 protein detects endogenous levels of MAPK3
Storage and Stability:	-20°C/1 year