

Antibodies Product Datasheet

Product Name:	HP-1 γ mouse Monoclonal Antibody(3B9)
Catalog Number:	N1685
Reactivity:	Human,Mouse,Rat
Applications:	WB,IHC-p,IF(paraffin section)
Dilutions:	WB 1:500-2000,IHC-p 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of HP-1 γ
Concentration:	1 mg/ml
Gene Name:	CBX3
Protein Name:	Chromobox protein homolog 3 (HECH) (Heterochromatin protein 1 homolog gamma) (HP1 gamma) (Modifier 2 protein)
Other Name:	Chromobox protein homolog 3 (HECH) (Heterochromatin protein 1 homolog gamma) (HP1 gamma) (Modifier 2 protein)
Human Swiss Prot No.:	Q13185
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous HP-1 γ protein
Storage and Stability:	-20°C/1 year