

Antibodies Product Datasheet

Product Name:	ATG5 mouse Monoclonal Antibody(3C7)
Catalog Number:	N1695
Reactivity:	Human,Mouse,Rat
Applications:	IF/ICC, WB, IHC-pIP
Dilutions:	IF: 1:50-200 IP: 1:50-200 WB 1:500-2000, IHC-p 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of ATG5
Concentration:	1 mg/ml
Gene Name:	ATG5 APG5L ASP
Protein Name:	Autophagy protein 5 (APG5-like) (Apoptosis-specific protein)
Other Name:	Autophagy protein 5 (APG5-like) (Apoptosis-specific protein)
Human Swiss Prot No.:	Q9H1Y0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous ATG5 protein
Storage and Stability:	-20°C/1 year