

Antibodies Product Datasheet

Product Name:	Acetyl NF kB P65(K314/K315) mouse Monoclonal Antibody(5G11)
Catalog Number:	N1698
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic Peptide of Acetyl NF kB P65(K314/K315)
Concentration:	1 mg/ml
Gene Name:	RELA NFKB3
Protein Name:	Transcription factor p65 (Nuclear factor NF-kappa-B p65 subunit) (Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3)
Other Name:	Transcription factor p65 (Nuclear factor NF-kappa-B p65 subunit) (Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3)
Human Swiss Prot No.:	Q04206
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous Acetyl NF kB P65(K314/K315) protein
Storage and Stability:	-20°C/1 year