

Antibodies Product Datasheet

Product Name:	Glu-Glu Tag mouse Monoclonal Antibody(2G1)
Catalog Number:	N1718
Reactivity:	Species independent
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Glu-Glu Tag
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	N/A
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects Glu-Glu Tag protein
Storage and Stability:	-20°C/1 year