

Antibodies Product Datasheet

Product Name:	Collagen II mouse Monoclonal Antibody(7F9)
Catalog Number:	N1719
Reactivity:	Human,Mouse,Rat
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic Peptide of Collagen II
Concentration:	1 mg/ml
Gene Name:	COL2A1
Protein Name:	Collagen alpha-1(II) chain (Alpha-1 type II collagen) [Cleaved into: Collagen alpha-1(II) chain; Chondrocalcin]
Other Name:	Collagen alpha-1(II) chain (Alpha-1 type II collagen) [Cleaved into: Collagen alpha-1(II) chain; Chondrocalcin]
Human Swiss Prot No.:	P02458
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous Collagen II protein
Storage and Stability:	-20°C/1 year