

Antibodies Product Datasheet

Product Name:	CD146 mouse mAb(PT0278)
Catalog Number:	N1730
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section),WB
Dilutions:	IHC-p 1:100-300, WB 1:500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CD146
Concentration:	1 mg/ml
Gene Name:	MCAM MUC18
Protein Name:	CD146
Other Name:	Cell surface glycoprotein MUC18 (Cell surface glycoprotein P1H12) (Melanoma cell adhesion molecule) (Melanoma-associated antigen A32) (Melanoma-associated antigen MUC18) (S-endo 1 endothelial-associated antigen) (CD antigen CD146)
Human Swiss Prot No.:	P43121
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CD146
Storage and Stability:	-20°C/1 year