

Antibodies Product Datasheet

Product Name:	ERCC1 mouse mAb(PT0756)
Catalog Number:	N1736
Reactivity:	Human
Applications:	WB,IF/ICC
Dilutions:	WB 1:500-2000, IF 1:100-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgG1, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ERCC1
Concentration:	1 mg/ml
Gene Name:	ERCC1
Protein Name:	ERCC1
Other Name:	DNA excision repair protein ERCC-1
Human Swiss Prot No.:	P07992
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ERCC1
Storage and Stability:	-20°C/1 year