

Antibodies Product Datasheet

Product Name:	Alpha-1-Fetoprotein(AFP) mouse mAb(PT0679)
Catalog Number:	N1737
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgG1, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Alpha-1-Fetoprotein(AFP)
Concentration:	1 mg/ml
Gene Name:	AFP
Protein Name:	Alpha-1-Fetoprotein(AFP)
Other Name:	Alpha-fetoprotein (Alpha-1-fetoprotein) (Alpha-fetoglobulin)
Human Swiss Prot No.:	P02771
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Alpha-1-Fetoprotein(AFP)
Storage and Stability:	-20°C/1 year