

Antibodies Product Datasheet

Product Name:	Postmeiotic Segregation Increased 2(PMS2) mouse mAb(PT2116)
Catalog Number:	N1752
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC-p 1:100-500, IF 1:100-500
Modification:	Unmodified/ Total
Source:	Mouse,Monoclonal
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Postmeiotic Segregation Increased 2(PMS2)
Concentration:	
Gene Name:	PMS2 PMSL2
Protein Name:	Postmeiotic Segregation Increased 2(PMS2)
Other Name:	Mismatch repair endonuclease PMS2 (EC 3.1.-.-) (DNA mismatch repair protein PMS2) (PMS1 protein homolog 2)
Human Swiss Prot No.:	P54278
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.151% sodium azide.
Specificity:	This antibody detects endogenous levels of Postmeiotic Segregation Increased 2(PMS2) at Human
Storage and Stability:	-20°C/1 year