

Antibodies Product Datasheet

Product Name:	Tissue Inhibitor of Metalloproteinases 1(PT2236) mouse mAb
Catalog Number:	N1865
Reactivity:	Human, Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	WB 500-2000 IHC-p 1:100-500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgG1, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Tissue Inhibitor of Metalloproteinases 1(TIMP1)
Concentration:	
Gene Name:	TIMP1 CLGI TIMP
Protein Name:	Metalloproteinase inhibitor 1 (Erythroid-potentiating activity) (EPA) (Fibroblast collagenase inhibitor) (Collagenase inhibitor) (Tissue inhibitor of metalloproteinases 1) (TIMP-1)
Other Name:	
Human Swiss Prot No.:	P01033
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Metalloproteinase inhibitor 1 (Erythroid-potentiating activity) (EPA) (Fibroblast collagenase inhibitor) (Collagenase inhibitor) (Tissue inhibitor of metalloproteinases 1) (TIMP-1). Heat-induced epitope ret
Storage and Stability:	-20°C/1 year