

Antibodies Product Datasheet

Product Name:	Napsin A mouse mAb(PT2252)
Catalog Number:	N1876
Reactivity:	Human, Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC-p 1:100-500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgG1, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Napsin A
Concentration:	
Gene Name:	NAPSA NAP1 NAPA
Protein Name:	Napsin A
Other Name:	Napsin-A (EC 3.4.23.-) (Aspartyl protease 4) (ASP4) (Asp 4) (Napsin-1) (TA01/TA02)
Human Swiss Prot No.:	O96009
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Napsin A. Heat-induced epitope retrieval (HIER) TRIS-EDTA of pH9.0 was highly recommended as antigen repair method in paraffin section
Storage and Stability:	-20°C/1 year