

## Antibodies Product Datasheet

|                               |                                                                                                                                                                                                             |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Neuron-Specific Enolase(NSE) mouse mAb(PT2281)                                                                                                                                                              |
| <b>Catalog Number:</b>        | N1881                                                                                                                                                                                                       |
| <b>Reactivity:</b>            | Human, Mouse                                                                                                                                                                                                |
| <b>Applications:</b>          | IHC-p,IF(paraffin section)                                                                                                                                                                                  |
| <b>Dilutions:</b>             | IHC-p 1:100-500                                                                                                                                                                                             |
| <b>Modification:</b>          | Unmodified/ Total                                                                                                                                                                                           |
| <b>Source:</b>                | Monoclonal Mouse/IgG1, Kappa                                                                                                                                                                                |
| <b>Purification:</b>          | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.                                                                                                  |
| <b>Immunogen:</b>             | Synthesized peptide derived from human Neuron-Specific Enolase(NSE)                                                                                                                                         |
| <b>Concentration:</b>         |                                                                                                                                                                                                             |
| <b>Gene Name:</b>             | ENO2                                                                                                                                                                                                        |
| <b>Protein Name:</b>          | Neuron-Specific Enolase(NSE)                                                                                                                                                                                |
| <b>Other Name:</b>            | Gamma-enolase (EC 4.2.1.11) (2-phospho-D-glycerate hydro-lyase) (Enolase 2) (Neural enolase) (Neuron-specific enolase) (NSE)                                                                                |
| <b>Human Swiss Prot No.:</b>  | P09104                                                                                                                                                                                                      |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                     |
| <b>Specificity:</b>           | This antibody detects endogenous levels of human Neuron-Specific Enolase(NSE). Heat-induced epitope retrieval (HIER) TRIS-EDTA of pH9.0 was highly recommended as antigen repair method in paraffin section |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                                                                |