

Antibodies Product Datasheet

Product Name:	β -tubulin (HRP Conjugated) Monoclonal Antibody(5G3)
Catalog Number:	N2000
Reactivity:	Human, Rat, Mouse, Mk, Dg, Ch, Hamster, Rabbit, sheep, Insect, Yeast
Applications:	WB, IF, IHC
Dilutions:	WB 1:1000-1:3000,IHC 1:50-1:200 IF 1:200
Modification:	Unmodified/Total
Source:	Monoclonal, Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of β -tubulin (HRP Conjugated)
Concentration:	1 mg/ml
Gene Name:	TUBB
Protein Name:	Tubulin β
Other Name:	Tubulin beta chain (Tubulin beta-5 chain)
Human Swiss Prot No.:	P07437
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The β -tubulin antibody can detects endogenous β -tubulin protein.
Storage and Stability:	Store at -20° C/1 year