

Antibodies Product Datasheet

Product Name:	MEK5 Mouse mAb(2G10)
Catalog Number:	N2005
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Mouse,monoclonal
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human MEK5
Concentration:	1 mg/ml
Gene Name:	MAP2K5 MEK5 MKK5 PRKMK5
Protein Name:	Kininogen 1
Other Name:	Dual specificity mitogen-activated protein kinase kinase 5 (MAP kinase kinase 5) (MAPKK 5) (EC 2.7.12.2) (MAPK/ERK kinase 5) (MEK 5)
Human Swiss Prot No.:	Q13163
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MEK5 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year