

## Antibodies Product Datasheet

<b>Product Name:</b>	ACE2 Mouse mAb(2C1)
<b>Catalog Number:</b>	N2006
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC, IF
<b>Dilutions:</b>	WB 1:500-2000 IHC 1:50-200. IF 1:50-200
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Mouse,monoclonal
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human ACE2
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	ACE2 UNQ868/PRO1885
<b>Protein Name:</b>	ACE2
<b>Other Name:</b>	Angiotensin-converting enzyme 2 (EC 3.4.17.23) (ACE-related carboxypeptidase) (Angiotensin-converting enzyme homolog) (ACEH) (Metalloprotease MPROT15) [Cleaved into: Processed angiotensin-converting enzyme 2]
<b>Human Swiss Prot No.:</b>	Q9BYF1
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of ACE2 at Human, Mouse,Rat
<b>Storage and Stability:</b>	Store at -20° C/1 year