

Antibodies Product Datasheet

Product Name:	Zip9 Mouse mAb(Mix-mA)
Catalog Number:	N2009
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Mouse,monoclonal
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Zip9
Concentration:	1 mg/ml
Gene Name:	SLC39A9 ZIP9 UNQ714/PRO1377
Protein Name:	ZIP9
Other Name:	Zinc transporter ZIP9 (Solute carrier family 39 member 9) (Zrt- and Irt-like protein 9) (ZIP-9)
Human Swiss Prot No.:	Q9NUM3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Zip9 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year