

Antibodies Product Datasheet

Product Name:	DDDDK-Tag(binds to flag sequence) mouse mAb(12M11)
Catalog Number:	N2010
Reactivity:	Species independent
Applications:	WB, ELISA, IP, IF, CoIP
Dilutions:	WB 1:10000-50000 ELISA 1:20000-50000 IP 1:50-1:200 IF 1:1000 -3000
Modification:	Unmodified/Total
Source:	Mouse VH IgG1 and VL Kappa was expressed from 293F cells
Purification:	Protein A purified
Immunogen:	Synthetic Peptide of Flag-Tag
Concentration:	
Gene Name:	Flag tag; Flag-tag,DDDDK TAG, DDDDK-TAG, DYKDDDDK tag,DYKDDDDK-tag
Protein Name:	Flag-tag
Other Name:	
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	This recombinant antibody can highly specifically recognize the Flag tags of C-terminal and N-terminal of recombinant protein. Variable region gene of this antibody comes from the mouse gene library immunized with DDDK peptide. Constant region of this ant
Storage and Stability:	Store at -20° C/1 year