

Antibodies Product Datasheet

Product Name:	CD206 (PTR2276) mouse mAb
Catalog Number:	N2013
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF, ELISA
Dilutions:	WB 1:500-2000. IF 1:100-500. ELISA 1:1000-5000
Modification:	Unmodified/Total
Source:	Mouse, Monoclonal/IgG1, kappa
Purification:	Protein G
Immunogen:	Synthesized peptide derived from human protein.AA range:750-850
Concentration:	
Gene Name:	MRC1
Protein Name:	CD206
Other Name:	MRC1; CLEC13D; CLEC13DL; MRC1L1; Macrophage mannose receptor 1; MMR; C-type lectin domain family 13 member D; C-type lectin domain family 13 member D-like; Macrophage mannose receptor 1-like protein 1; CD206
Human Swiss Prot No.:	P22897
Formulation:	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Specificity:	This antibody detects endogenous levels of CD206.
Storage and Stability:	Store at -20° C/1 year