

Antibodies Product Datasheet

Product Name:	ERdj3 (PTR1359) mouse mAb
Catalog Number:	N2015
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF, ELISA
Dilutions:	WB 1:500-2000. IF 1:100-500. ELISA 1:1000-5000
Modification:	Unmodified/Total
Source:	Mouse, Monoclonal/IgG
Purification:	Protein G
Immunogen:	Synthesized peptide derived from human protein.AA range: 1-100
Concentration:	
Gene Name:	DNAJB11
Protein Name:	ERdj3
Other Name:	DNAJB11; EDJ; ERJ3; HDJ9; PSEC0121; DnaJ homolog subfamily B member 11; APOBEC1-binding protein 2; ABBP-2; DnaJ protein homolog 9; ER-associated DNAJ; ER-associated Hsp40 co-chaperone; ER-associated dnaJ protein 3; ERdj3; ERj3p; HEDJ; Human Q9UBS4
Human Swiss Prot No.:	
Formulation:	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Specificity:	This antibody detects endogenous levels of ERdj3.
Storage and Stability:	Store at -20° C/1 year