

Antibodies Product Datasheet

Product Name:	Perforin (ABT365) mouse mAb
Catalog Number:	N2016
Reactivity:	Human
Applications:	IHC, IF, ELISA
Dilutions:	IHC 1:50-200. IF 1:50-200. ELISA 1:500-5000
Modification:	Unmodified/Total
Source:	Mouse, Monoclonal/IgG1, kappa
Purification:	The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Perforin AA range: 1-100
Concentration:	
Gene Name:	PRF1 PFP
Protein Name:	Perforin
Other Name:	Cytolysin; FLH2; HPLH2; Lymphocyte pore-forming protein; P1; PERF_HUMAN; perforin 1 (pore forming protein); Perforin 1; Perforin-1; PFP; PGFL; PIGF; PIGF-2; PLGF; Pore forming protein; prf1; SHGC-10760
Human Swiss Prot No.:	P14222
Formulation:	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Specificity:	The antibody can specifically recognize human Perforin protein.
Storage and Stability:	Store at -20° C/1 year