

Antibodies Product Datasheet

Product Name:	MutL Protein Homolog 1(MLH1) (ABT-MLH1) mouse mAb
Catalog Number:	N2017
Reactivity:	Human
Applications:	IHC, ELISA
Dilutions:	IHC 1:50-100. ELISA 1:500-5000
Modification:	Unmodified/Total
Source:	Mouse, Monoclonal/IgG1, kappa
Purification:	The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MutL Protein Homolog 1(MLH1) AA range: 400-500
Concentration:	
Gene Name:	MLH1 COCA2
Protein Name:	MLH1
Other Name:	
Human Swiss Prot No.:	P40692
Formulation:	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Specificity:	This antibody detects endogenous levels of MLH1 protein.
Storage and Stability:	Store at -20° C/1 year