

Antibodies Product Datasheet

Product Name:	CD74 (ABT202) mouse mAb
Catalog Number:	N2018
Reactivity:	Human
Applications:	IHC, ELISA
Dilutions:	IHC 1:200-400. ELISA 1:500-5000
Modification:	Unmodified/Total
Source:	Mouse, Monoclonal/IgG2a, kappa
Purification:	The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CD74 AA range: 150-250
Concentration:	
Gene Name:	CD74 DHLAG
Protein Name:	CD74
Other Name:	CD 74 antibody; CD74 antibody; CD74 antigen (invariant polypeptide of major histocompatibility complex; class II antigen-associated) antibody; CD74 antigen antibody; CD74 molecule antibody; CD74 molecule; major histocompatibility complex; class II invaria
Human Swiss Prot No.:	P04233
Formulation:	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Specificity:	This antibody detects endogenous levels of CD74 protein.
Storage and Stability:	Store at -20° C/1 year