

Antibodies Product Datasheet

Product Name:	MutS Protein Homolog 2(MSH2) (ABT-MSH2) mouse mAb
Catalog Number:	N2020
Reactivity:	Human
Applications:	IHC, ELISA
Dilutions:	IHC 1:50-100. ELISA 1:500-5000
Modification:	Unmodified/Total
Source:	Mouse, Monoclonal/IgG1, kappa
Purification:	The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MutS Protein Homolog 2(MSH2) AA range: 600-700
Concentration:	
Gene Name:	MSH2
Protein Name:	MSH2
Other Name:	DNA mismatch repair protein Msh2 (hMSH2; MutS protein homolog 2)
Human Swiss Prot No.:	P43246
Formulation:	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Specificity:	This antibody detects endogenous levels of MSH2 protein.
Storage and Stability:	Store at -20° C/1 year