

Antibodies Product Datasheet

Product Name:	p120 catenin (ABT-p120) mouse mAb
Catalog Number:	N2021
Reactivity:	Human, Mouse, Rat
Applications:	IHC, ELISA
Dilutions:	IHC 1:200-400. ELISA 1:500-5000
Modification:	Unmodified/Total
Source:	Mouse, Monoclonal/IgG2b, kappa
Purification:	The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human p120 catenin AA range: 600-700
Concentration:	
Gene Name:	CTNND1 KIAA0384
Protein Name:	p120 Catenin
Other Name:	Catenin delta-1 (Cadherin-associated Src substrate; CAS; p120 catenin; p120(ctn); p120(cas))
Human Swiss Prot No.:	O60716
Formulation:	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Specificity:	The antibody can specifically recognize human p120 catenin protein.
Storage and Stability:	Store at -20° C/1 year