

Antibodies Product Datasheet

Product Name:	Ki-67 (ABT001) mouse mAb
Catalog Number:	N2022
Reactivity:	Human
Applications:	IHC, ELISA
Dilutions:	IHC 1:200-400. ELISA 1:500-5000
Modification:	Unmodified/Total
Source:	Mouse, Monoclonal/IgG2b, kappa
Purification:	The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Ki-67 AA range: 1200-1400
Concentration:	
Gene Name:	MKI67
Protein Name:	Ki-67
Other Name:	MKI67; Antigen KI-67
Human Swiss Prot No.:	P46013
Formulation:	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Specificity:	The antibody can specifically recognize human Ki-67 protein.
Storage and Stability:	Store at -20° C/1 year