

Antibodies Product Datasheet

Product Name:	Lamin A/C (phospho Ser392) (PTR1133) mouse mAb
Catalog Number:	N2025
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IF, ELISA
Dilutions:	WB 1:500-2000. IF 1:100-500. ELISA 1:1000-5000
Modification:	Unmodified/Total
Source:	Mouse, Monoclonal/IgG1, kappa
Purification:	Protein G
Immunogen:	Synthesized peptide derived from human protein. AA range: 350-450
Concentration:	1 mg/ml
Gene Name:	LMNA LMN1
Protein Name:	Lamin A/C (phospho Ser392)
Other Name:	Prelamin-A/C [Cleaved into: Lamin-A/C (70 kDa lamin) (Renal carcinoma antigen NY-REN-32)]
Human Swiss Prot No.:	P02545
Formulation:	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Specificity:	This antibody detects endogenous levels of Lamin A/C (phospho Ser392) protein.
Storage and Stability:	Store at -20° C/1 year