

Antibodies Product Datasheet

Product Name:	eIF2 α (Phospho Ser51) (PTR1316) mouse mAb
Catalog Number:	N2026
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF, ELISA
Dilutions:	WB 1:500-2000. IF 1:100-500. ELISA 1:1000-5000
Modification:	Unmodified/Total
Source:	Mouse, Monoclonal/IgG2b, kappa
Purification:	Protein G
Immunogen:	Synthesized peptide derived from human eIF2 α AA range: 1-100
Concentration:	1 mg/ml
Gene Name:	EIF2S1 EIF2A
Protein Name:	EIF2S1
Other Name:	
Human Swiss Prot No.:	P05198
Formulation:	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Specificity:	This antibody detects endogenous levels of eIF2 α (Phospho Ser51) protein.
Storage and Stability:	Store at -20° C/1 year